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## Pedikelschraubensystem

### PsOne™/PsTwo™

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## PsOne™/PsTwo™ pedicle screw system

### Intended Purpose

The ARCA-MEDICA GmbH long term implantable PsOne™/PsTwo™ pedicle screw system is comprised of a range of pedicle screws and rods that are used for correction of spinal deformities, fractures, degenerative and dislocated vertebrae in the thoracic, sacral and lumbar region.

All implants are only intended for single use and must not be reused under any circumstances.

### Clinical benefit

- Correction of spinal deformities, fractures, degenerative and dislocated vertebrae in the thoracic, sacral and lumbar region
- Stabilization of vertebrae and segments of the vertebral column in the thoracic, sacral and lumbar region
- pain relief and improvement of the quality of life

### INDICATIONS

- Degenerative disc disease (DDD) (defined as discogenic back pain with degeneration of the intervertebral disc confirmed by medical history and imaging examination)
- Spondylolisthesis
- Trauma (i.e. fracture or dislocation)
- Spinal stenosis
- Spondylitis
- Significant mechanical instability or deformity of the thoracic, lumbar and sacral spine that leads to secondary severe spondylolisthesis of the L5-S1 vertebrae

### CONTRAINDICATIONS

- Contraindications include, but are not limited to:
- Absent osseous structures that make sound anchoring of the implant impossible (e.g. osteoporosis or bone decalcification)
  - Local infection in the immediate vicinity of the implant site
  - Inadequately controlled systemic infections
  - Fever or leukocytosis
  - Obesity
  - Pregnancy
  - Mental illnesses
  - Severely distorted anatomy due to congenital anomalies
  - Medical or surgical conditions that prevent the possible benefits of surgery on the spinal column
  - Suspected or documented allergy or intolerance to the implant material
  - Any case in which the use of a pedicle screw system is not necessary
  - Any case in which the implant components selected for use are too big or small to achieve a successful outcome
  - Any patient with inadequate tissue cover over the surgical site or inadequate bone mass or quality at the surgical site
  - Any patient in which the implant would impair the anatomical structures or the expected physiological performance
  - Paediatric patients or in case of general skeletal growth
  - A lack of willingness of the patient to follow the post-surgical instructions
  - Any case that is not described in the indications

- All set screws must be tightened as described in the surgical technique of the PsOne™/PsTwo™ pedicle screw system.

### POSTOPERATIVE MANAGEMENT

In the postoperative phase, pay particular attention to the supply of individual patient information by the treating physician.

- The aftercare is defined by the treating physician and is individually tailored to the patient. The patient must be given detailed instructions on the use and restrictions of the system
- The patient must limit their physical activity during the healing process. This includes, in particular, the lifting of loads, rotating movements and any type of sports activity. Extreme stresses, falls, vibrations, sudden shocks or impacts on the spinal column should be avoided.
- The patient must be informed that bending, loosening and/or fracture of the system are complications that can occur as a result of excessive or premature stress or muscle activity.
- During the healing process, it is recommended that the patient does not smoke, does not use any nicotine products, does not drink alcohol and does not take any non-steroidal or anti-inflammatory medication, such as aspirin.
- Ensure that the patient is aware of the permanent physical limitation caused by spinal fusion and learns to compensate for this.
- It is important that the surgical site of the spinal column is kept immobilised until permanent bone fusion is established and confirmed by a radiological examination.
- An implant failure may necessitate revision surgery.
- The patient must attend regular medical follow-up monitoring of the implant.
- Removed implants are considered explants and must be handled in a manner that prevents their reuse in another surgical intervention. The components of the PsOne™/PsTwo™ pedicle screw system must never be reused under any circumstances.

### REVISION

The PsOne™/PsTwo™ pedicle screw system may only be removed with the specified Ps instruments. The explant procedure is described in the surgical technique.

### PACKAGING

 Sterilised by radiation

 Do not re-sterilise

The PsOne™/PsTwo™ implants are delivered sterile and ready for use. The content is sterile, unless the packaging is damaged, open or the expiration date on the label has expired. All instruments are delivered non-sterile.

 Do not reuse

All implants are only intended for single use and must not be reused under any circumstances. The reuse of a disposable device that has been in contact with blood, bone, tissue, other bodily fluids or other biological non-patient-related burden can injure the patient or user. Possible risks associated with the reuse of a disposable device include, among other things,

- Medication or drug abuse or alcohol addiction

Although not considered absolute contraindications, the following diseases must be considered potential factors for the non-use of this system:

- Severe bone resorption
- Osteomalacia

### POSSIBLE COMPLICATIONS

- Wear, bending or fracture of the implant components
- Loosening of one or all components
- Loss of fixation, dislocation, sintering
- Bursitis
- Foreign body reaction or allergic reaction to implants
- Post-operative change to the spinal curvature, loss of correction, height and/or reduction
- Infection
- Pseudomeningocele, fistula, persistent CSF leakage, meningitis
- Loss of neurological function, including paralysis and/or the development or persistence of pain, spasms, loss of sensation, sensory disorders such as tingling and/or impaired vision
- Arachnoiditis and/or muscular atrophy
- Urinary retention or loss of bladder control or other types of impairments of the urinary tract system
- Fracture, microfracture, resorption, damage or penetration of the bones of the spinal column (including the sacrum, pedicle and/or vertebral bodies)
- Soft tissue injuries, vascular or visceral injuries
- Nucleus pulposus herniation, disc herniation or degeneration at, above or below the surgery level
- Pseudoarthrosis, delayed fusion, malformation
- Inhibition of the potential growth of the operated spinal column section
- Loss of mobility or function of the spinal column
- Bone loss or reduced bone density, potentially caused by stress shielding
- Ileus, gastritis, constipation, loss of bowel control or other restrictions of the gastrointestinal system
- Bleeding, haematomas, occlusion, seroma, oedemas, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to the blood vessels or other impairments of the cardiovascular system
- Impairment of the reproductive system, including sterility, loss of consortium and sexual dysfunction
- Inability to carry out everyday activities
- Change in mental condition
- Death

### Patient group

The patient group for the PsOne™/PsTwo™ pedicle screw system includes skeletally mature patients requiring immobilization and stabilization of spinal segments through posterior, non-cervical pedicular fixation of the thoracic, lumbar, and sacral spine.

### Intended User

The PsOne™/PsTwo™ pedicle screw system is intended to be implanted by adequately trained surgeon for spinal surgery and who is familiar with the biomechanical conditions of the spinal

mechanical failure, serious injuries, transmission of infectious diseases and death.

 Damaged packaging

The packaging of the individual components must be intact upon receipt. All packaging must be carefully checked for damage before use. Damaged packaging or devices must not be used and must be returned.

 Expires on

Always check the expiry date of the device prior to use. Expired devices must not be used.

All instruments must be cleaned, disinfected and sterilised as specified in the reprocessing instructions.

### STORAGE CONDITIONS

 Always store the implant in the original protective packaging. Always remove the implant from the protective packaging only directly before use. Store the implant in a dry and dust-free place (standard hospital conditions). Avoid direct sunlight.

### MRI SAFETY

#### MR conditional

The components of the PsOne™/PsTwo™ pedicle screw system have not been tested for safety and compatibility in MR environments (magnetic resonance). They have not been tested for heating, migration or image artefacts in MR environments. The material, design and intended use of the PsOne™/PsTwo™ System is not suspected to induce any critical negative effects when subjecting patients an MRI treatment with 1.5 or 3 Tesla. However, the risk assessment must always be carried out by the physician on case-by-case basis. ARCA-MEDICA GmbH accepts no responsibility for the safety of the PsOne™/PsTwo™ pedicle screw system in MR environments.

### PRODUCT COMPLAINTS

ARCA-MEDICA GmbH must immediately be notified of any apparent (possible) malfunction. The article number(s), batch number(s) and a contact address must be provided with every complaint. ARCA-MEDICA GmbH must immediately be informed in the event of a (possible) malfunction that leads to/contributes to serious injury or the death of the patient. [www.arca-medica.com/unternehmen/kontakt/](http://www.arca-medica.com/unternehmen/kontakt/)

Driescheider Weg 40-42  
D-57610 Altenkirchen  
Germany  
Tel: +49 (0) 7223 80805 210  
Fax: +49 (0) 7223 80805 211  
[www.arca-medica.com](http://www.arca-medica.com)

column and the relevant surgical techniques. The use of these products must be in a clinical facility only.

### SAFETY INSTRUCTIONS

- Read these instructions carefully before using the device and familiarise yourself with the surgical technique.
- Store the instructions for use and surgical technique in an accessible place for all personnel.
- The surgeon must have theoretical as well as practical command of the recognised surgical techniques and is responsible for proper implantation.
- The manufacturer accepts no responsibility for complications due to incorrect diagnosis, incorrect implant selection, incorrect surgical techniques, limits of the treatment method or inadequate asepsis.
- Damage to the force-transmitting structures can lead to implant fracture and/or implant loss and cause dislocation and migration as well as other complications. To detect factors that facilitate any malfunction of the implant at an early stage, the implant should be regularly monitored after surgery using appropriate methods. Radiological monitoring of the correct position and integrity of the implant should take place in case of the sudden onset of severe pain or neurological symptoms.

### WARNINGS

- Do not use damaged or defective instruments! Request a repair or replacement from the manufacturer or distributor.
- The improper use of an instrument may damage/destroy the instrument and injure the patient and/or surgical staff.
- Only use original ARCA-MEDICA instruments. The use of other instruments may impede the successful outcome of the surgery and lead to injuries and/or failure of the implant.
- Never use implants if the packaging is damaged. An implant taken from damaged packaging is not sterile and/or is damaged and must therefore not be used.
- Never use implants whose expiration date has passed.
- The PsOne™/PsTwo™ pedicle screw system must not be used with components of other systems or manufacturers.
- All implants are only intended for single use and must not be reused under any circumstances. Do not reprocess and/or sterilise.

### IMPLANT SELECTION

Selection of the correct size of the PsOne™/PsTwo™ pedicle screw for every patient and vertebra is critical for the success of the surgery. Incorrectly dimensioned screws lead to complications, such as nerve damage, bleeding or implant failure.

### MATERIAL

The PsOne™/PsTwo™ pedicle screw system is made from an implant grade titanium alloy (Ti6Al4V ELI).

### INSTRUMENTS

The associated Ps instruments are made from various stainless steels, titanium alloys and polymers. All materials used are biocompatible.

### WARRANTY

Any product warranty expires:

- In case of damage caused by improper storage, handling, cleaning and/or sterilisation.
- The manufacturer accepts no responsibility for complications due to incorrect diagnosis, incorrect implant selection, incorrect surgical techniques, limits of the treatment method or inadequate asepsis.

### TRACEABILITY

The self-adhesive identification labels contained in the packaging must be stored to ensure traceability and enclosed with the case report and patient file.

The hospital must keep a list of devices in storage, including the article number and batch number, so that measures can quickly be taken in case of changes or recalls.

When passing on an ARCA-MEDICA product (with or without charge), it must be ensured that traceability (LOT tracking) is possible at all times.

### DISPOSAL


Implants that are no longer usable can be disposed of in accordance with local regulations or returned to the manufacturer free of charge for proper disposal. ARCA-MEDICA implants may only be returned if they are in their undamaged original packaging or have been cleaned, disinfected and sterilised. Proof of this must be attached to the outer packaging.

### ADDITIONAL INFORMATION

If you require or desire any additional information, please contact ARCA-MEDICA GmbH.

### SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of safety and clinical performance is available on request. Please visit [www.arca-medica.com/unternehmen/kontakt/](http://www.arca-medica.com/unternehmen/kontakt/).

	Caution	Refers to the need for the user to review the instructions for use for important safety information, such as warnings and precautions, which cannot be attached to the medical device itself for various reasons.
	Keep away from sunlight	Refers to a medical device that must be kept away from light sources.
	Store in a dry place	Refers to a medical device that must be protected against moisture.
	MR conditional	An object that is considered safe in an MR environment under certain conditions

### COMBINATION AND COMPATIBLE PRODUCTS

The stand-alone usage of the PsOne™/PsTwo™ pedicle screw system is only allowed for single level application. Starting with two levels and more, the usage of a cross connector is mandatory.

The PsOne™/PsTwo™ pedicle screw system is also intended to be combined with cages to further support the stabilisation.

The list of allowed products is part of the surgical technique document.

### PREOPERATIVE MANAGEMENT

- Only patients who satisfy the criteria described in the "Indications" section are selected.
- Patient information on the advantages and disadvantages of an implant and alternative treatment options as well as on possible complications and risks of the surgery.
- Handle and store the implant components with care. The implants must not be scratched or otherwise damaged. Implants and instruments must be protected during storage.
- At the time of the surgery, an adequate stock of implants in excess of the expected quantity must be available.
- The surgeon must have theoretical as well as practical command of the recognised surgical techniques and is responsible for proper implantation.
- The surgeon must be an adequately trained specialist for spinal surgery and be familiar with the biomechanical conditions of the spinal column and the relevant surgical techniques.
- The surgeon and surgical team must be intimately familiar with the surgical technique as well as with the intended implants and instruments.
- The surgeon must personally check whether the necessary components are available prior to surgery before using the system. The PsOne™/PsTwo™ pedicle screw system must not be combined with components other those specified in the surgical technique
- All instruments must be cleaned and sterilised before use.
- All instruments must be completely available and functional. Damaged or defective instruments must not be used and a replacement must be requested from the distributor or ARCA-MEDICA.

### INTRAOPERATIVE MANAGEMENT

- Extreme caution is required in the area of the spinal cord and nerve roots. Nerve damage can lead to the loss of neurological functions.
- The fracture, slipping or improper use of instruments or implant components can injure the patient or surgical staff.
- In particular, make sure that the implant surfaces are not scratched or scored, as such damage impairs the functional strength of the construct. Only use intact implants.
- Use an imaging technique to check the correct position of the system.
- Avoid the use of excess force when inserting the implant in order to protect the adjacent vertebral bodies and nerves.
- The rods must not be bent back and forth several times at the same point.

### DESCRIPTION OF THE SYMBOLS

Symbol	Symbol title	Symbol description
	Batch code	Indicates the manufacturer's batch designation so that the batch or lot can be identified.
	Article number	Shows the manufacturer's order number so that the medical device can be identified.
	Sterilised by radiation	Indicates a medical device that has been sterilised by radiation.
	Do not reuse	Refers to a medical device that is intended for single use or use on a single patient for a single treatment.
	Do not re-sterilise	Indicates a medical device that must not be re-sterilised.
	Note the instructions for use	Refers to the need for the user to consult the instructions for use.
	Expires on	Indicates the date after which the medical device must no longer be used.
	Date of manufacture	Indicates the date on which the medical device was manufactured.
	Manufacturer	Indicates the manufacturer of the medical device as per Medical Device Regulation (EU) 2017/745
	Do not use if packaging is damaged	Indicates a medical device that should not be used if the packaging is damaged or open.
	CE marking	CE marking as per Medical Device Regulation (EU) 2017/745
	Medical device	Indicates that the relevant object is a medical device.
	Unique device identifier	Indicates an item that contains information on a unique device identifier (UDI).
	Double sterile barrier system	Indicates two sterile barrier systems.